


# DIAGNOSTIC INSTRUMENTS





**Duplex Stethoscope 01-001**  
super lightweight



**Duplex-Baby Stethoscope 01-002**  
super lightweight



**Anaestophon Stethoscope 01-003**  
Super Lightweight



**Clar Mead Mirror komplett 6 Volt 01-004**



**01-007**  
Ø 50mm



**Pinard Wooden Stethoscope 01-008** 18 cm



**Pinard Aluminium Stethoscope 01-009** 15 cm



**Weber Aesthesiometer 01-010** 9 cm



**Ziegler 01-005**  
Ø 90mm with head band without protection cap



**Ziegler 01-006**  
Ø 90mm with head band with protection cap



**01-011** 16cm  
**01-012** 20cm  
**01-013** 25cm

**Magill Catheter Intubating Forceps**



**Magill Catheter Intubating Forcep Gold**  
**01-014** 6.5" (2.5" shaft) , Infant , hole ( dia 2.00 mm)  
**01-015** 7.5" (3" shaft) , Infant , hole ( dia 2.30 mm)  
**01-016** 7.5" (4" shaft) , Infant , hole ( dia 2.60 mm)

**Magill Catheter Intubating Forcep GOLD**  
**01-017** 8" (4" shaft) , Child , hole ( dia 3.00 mm)  
**01-018** 8" (4" shaft) , Child , oval tip ( 5.38 mm X 4.15 mm)  
**01-019** 9" (5" shaft) , Adult , oval ( 6.80 mm x 5.00 mm)  
**01-020** 10" (5" shaft) , Adult , oval ( 9.60 mm x 7.70 mm)



Perkussionshämmer | Martillos de percusión  
 Marteaux-percuteurs | Martelli percussori



**Neuro-Duoflex Hammer**  
01-021 22cm



**Dejerine Hammer**  
01-022 20cm



**Dejerine Hammer**  
01-023 20cm



**Troemner Hammer**  
01-024 145 g  
01-025 175 g  
01-026 235 g  
24.5 cm  
spare rubber



**Troemner-Simplex Hammer**  
01-027 21cm



**Buck Hammer**  
01-028 18.5cm



**Buck Hammer**  
01-029 18cm



**Taylor Hammer**  
01-030 18cm  
01-031 20cm



**Taylor Hammer**  
01-032 18cm



**Babinsky Hammer**  
01-033 22 cm



**Babinsky Hammer**  
01-034 22 cm



**Rabiner Hammer**  
01-035 25.5 cm  
with needle & brush



**Berliner Hammer**  
01-036 20 cm




**Traube Hammer**  
01-037 16 cm



**Wartenberg Pinwheel**  
01-038 19 cm

**MC-INTOSH FIBER OPTIC LARYNGOSCOPE BLADE**




BLADE				
Size	A	B	C	Working Length
0	77mm	53mm	14mm	55mm
1	95mm	72mm	43mm	70mm
2	111mm	86mm	44mm	78mm
3	133mm	109mm	44mm	108mm
4	151mm	129mm	44mm	133mm

**Mc-Intosh Laryngoscope Fiber Optic**  
Halogen illumination. Manufactured from AISI 304 Stainless-Steel. Halogen illumination 2.5 V. The Bulb is inside the handle in order to simplify the sterilization of the blades. All blades and handles are interchangeable with other leading manufacturer. All sets can be combined as per your request, and delivered in Green plastic case.

**Blades Availability:**

- ★ Matt Finish
- ★ Mirror Finish
- ★ Sand Blast Finish
- ★ Pls, confirm in your order list

**MC-INTOSH LARYNGOSCOPE**



BLADE				
Size	A	B	C	Working Length
0	77mm	53mm	14mm	55mm
1	95mm	72mm	43mm	70mm
2	111mm	86mm	44mm	78mm
3	133mm	109mm	44mm	108mm
4	151mm	129mm	44mm	133mm

**Fiber Optic Laryngoscope**  
Plastic case, Halogen illumination. Manufactured from AISI 304 Stainless Steel Halogen Illumination 2.5 V. The Bulb is inside the handle in order to simplify the sterilization of the blades. All blades and handles are interchangeable with other leading manufacturer. All sets can be combined as per your request, and delivered in Green plastic case.

**Blades Availability:**

- ★ Matt Finish
- ★ Mirror Finish
- ★ Sand Blast Finish
- ★ Pls, confirm in your order list

**Mc-INTOSH FIBER OPTIC LARYNGOSCOPE SET HALOGEN ILLUMINATION**

- 01-039 PEDIATRIC SET**  
Set of 0,1,2 Blades  
Supplied in Green plastic case
- 01-040 STANDARD SET**  
Set of 1,2,3 Blades  
Supplied in Green plastic case
- 01-041 ADULT SET**  
Set of 1,2,3,4 Blades  
Supplied in Green plastic case



Laryngoscope Battery  
Handle for Adult  
C type 2 batteries  
**01-047**



Mc-Intosh Fiber Optic  
Laryngoscope Set  
Halogen Illumination  
Manufactured by AISI:  
304 Stainless Steel Material



Fig. 0  
**01-042**



Fig. 1  
**01-043**



Fig. 2  
**01-044**



Fig. 3  
**01-045**



Fig. 4  
**01-046**

**Mc-INTOSH FIBER OPTIC LARYNGOSCOPE SET HALOGEN ILLUMINATION**

- AVAILABLE SETS WITH BATTERY HANDLE**
- 01-048 PEDIATRIC SET**  
Set of 0,1,2 Blades  
Supplied in Green plastic case
  - 01-049 STANDARD SET**  
Set of 1,2,3 Blades  
Supplied in Green plastic case
  - 01-050 ADULT SET**  
Set of 1,2,3,4 Blades  
Supplied in Green plastic case



**01-056**  
Laryngoscope Battery  
Handle for Adult  
C type 2 Batteries



Fig. 0  
**01-051**



Fig. 1  
**01-052**



Fig. 2  
**01-053**



Fig. 3  
**01-054**



Fig. 4  
**01-055**



**01-057**  
Laryngoscope  
Bulb Small




**01-058**  
Laryngoscope  
Bulb Large

**Mc-Intosh Laryngoscope Set**  
of Fiber Optic  
Halogen Illumination  
Manufactured From: AISI 304  
Stainless Steel Material

Laryngoskope | Laringoscopios  
Laryngoscopes | Laringoscoopi

### FIBER OPTIC LARYNGOSCOPE MILLER



#### BLADE

Size	A	B	C	Working Length
0	75mm	53mm	26mm	55mm
1	111mm	80mm	26mm	80mm
2	151mm	134mm	28mm	130mm
3	188mm	167mm	28mm	170mm
4	203mm	180mm	28mm	180mm



**01-067**  
Laryngoscope Battery Handle for Adult  
C type 2 Batteries

### FIBER OPTIC LARYNGOSCOPE MILLER

#### AVAILABLE SETS WITH BATTERY HANDLE


**01-059 PEDIATRIC SET**  
Set of 0,1,2  
Blades  
Supplied in Green plastic case

**01-060 STANDARD SET**  
Set of 1,2,3  
Blades  
Supplied in Green plastic case

**01-061 ADULT SET**  
Set of 1,2,3,4  
Blades  
Supplied in Green plastic case



### MC-INTOSH LARYNGOSCOPE



#### BLADE

Size	A	B	C	Working Length
0	80mm	58mm	30mm	55mm
1	97mm	76mm	36mm	70mm
2	111mm	88mm	36mm	78mm
3	133mm	112mm	36mm	108mm
4	154mm	131mm	38mm	133mm

#### Mc-Intosh Laryngoscope

Plastic case, Standard illumination. Manufactured from AISI 304 Stainless Steel Matt Type, Standard Illumination 2.5 V. All blades and handles are interchangeable with other leading manufacturer. All sets can be combined as per your request, and delivered in Gray plastic case.

#### Blades Availability:

- \* Matt Finish
- \* Mirror Finish
- \* Sand Blast Finish
- \* Pls, confirm in your order list



**01-077**  
With one Spare Bulb



**01-076**  
Laryngoscope Battery Handle for Adult  
C type 2 Batteries



Mc-Intosh Laryngoscope Set  
Standard Illumination  
Manufactured From  
Stainless Steel Material

### CONVENTIONAL LARYNGOSCOPE SET

#### Mc-INTOSH LARYNGOSCOPE SET STANDARD ILLUMINATION

#### AVAILABLE SETS WITH BATTERY HANDLE


**01-068 PEDIATRIC SET**  
Set of 0,1,2 Blades  
Supplied in Gray plastic case

**01-069 STANDARD SET**  
Set of 1,2,3 Blades  
Supplied in Gray plastic case

**01-070 ADULT SET**  
Set of 1,2,3,4 Blades  
Supplied in Gray plastic case



**MC-INTOSH LARYNGOSCOPE**



**BLADE**

Size	A	B	C	Working Length
0	80mm	58mm	30mm	55mm
1	97mm	76mm	36mm	70mm
2	111mm	88mm	36mm	78mm
3	133mm	112mm	36mm	108mm
4	154mm	131mm	38mm	133mm

**Mc-Intosh Laryngoscope**

Plastic case, Standard illumination. Manufactured from AISI 304 Stainless Steel. All blades and handles are interchangeable with other leading manufacturer. All sets can be combined as per your request, and delivered in Blue plastic case.

**Blades Availability:**

- ★ Matt Finish
- ★ Mirror Finish
- ★ Sand Blast Finish
- ★ Pls, confirm in your order list



**01-087**  
With two Spare Bulbs



**01-086**  
Laryngoscope Battery Handle for Adult C type 2 batteries



**CONVENTIONAL LARYNGOSCOPE SET**

**Mc-INTOSH LARYNGOSCOPE CONVENTIONAL STANDARD ILLUMINATION**


**AVAILABLE SETS WITH BATTERY HANDLE**

- 01-078 PEDIATRIC SET**  
Set of 0,1,2 Blades.  
Supplied in Blue plastic case
- 01-079 STANDARD SET**  
Set of 1,2,3 Blades.  
Supplied in Blue plastic case
- 01-080 ADULT SET**  
Set of 1,2,3,4 Blades.  
Supplied in Blue plastic case



**Mc-Intosh Laryngoscope Conventional Standard Illumination**  
Manufactured From Stainless Steel material

**LARYNGOSCOPE MILLER**



**BLADE**

Size	A	B	C	Working Length
0	75mm	53mm	26mm	55mm
1	111mm	80mm	26mm	80mm
2	151mm	134mm	28mm	130mm
3	188mm	167mm	28mm	170mm
4	203mm	180mm	28mm	180mm

**Mc-Intosh Laryngoscope**

Plastic case, Standard illumination. Manufactured from AISI 304 Stainless Steel. All blades and handles are interchangeable with other leading manufacturer. All sets can be combined as per your request, and delivered in Blue plastic case.

**Blades Availability:**

- ★ Matt Finish
- ★ Mirror Finish
- ★ Sand Blast Finish
- ★ Pls, confirm in your order list



**01-097**  
With two Spare Bulbs



**01-096**  
Laryngoscope Battery Handle for Adult C type 2 batteries

**CONVENTIONAL LARYNGOSCOPE SET**

**LARYNGOSCOPE SET OF CONVENTIONAL MILLER STANDARD ILLUMINATION**

**AVAILABLE SETS WITH BATTERY HANDLE**

- 01-088 PEDIATRIC SET**  
Set of 0,1,2 Blades  
Supplied in Blue plastic case
- 01-089 STANDARD SET**  
Set of 1,2,3 Blades  
Supplied in Blue plastic case
- 01-090 ADULT SET**  
Set of 1,2,3,4 Blades  
Supplied in Blue plastic case



**Laryngoscope Set of conventional Miller Standard Illumination**  
Manufactured From Stainless Steel Material



Laryngoskope | Laringoscopios  
Laryngoscopes | Laringoscopi

**Laryngoscope Flexible (Mc-Coy)**



Available blade sizes 2,3 & 4 with 'AA' & 'C' type battery handles  
**01-100**  
**01-101**  
**01-102**

**01-099**  
'C' type handle

**01-098**  
**CONVENTIONAL SET**  
**Flexible Laryngoscope**  
**Conventional**  
Halogen Illumination  
Manufactured From  
Stainless Steel High  
grade 304 material.



Available Blade sizes 2,3 & 4 with 'AA' & 'C' type battery handles  
**01-105**  
**01-106**  
**01-107**

**01-103 FIBER OPTIC SET**  
**Flexible Laryngoscope**  
**Set of Fiber Optic**  
Halogen Illumination Manufactured  
From Stainless Steel High grade 304  
material packed in green plastic box.

Laryngoscope Battery  
Handle for Adult  
AA Batteries  
**01-104**



**DIAGNOSTIC PENTOSCOPE SET**  
**01-108**  
2.5 V Standard illumination for use with 2 batteries type AA.  
Consisting:  
Universal Handle thread type fit for all Pentoscope Heads like.

- Otoscope Head.
- Ophthalmoscope Head.
- Throat Branch Head.
- Mouth Mirrors Head.

It have two exploring mirrors No. 4 & 5 mm.  
Its spare 06 disposable Plastic ear Canulas are with sizes 2 x 2.50 mm, 2 x 3.50 mm & 2 x 4.5 mm.  
On special order it can be as 2 x 1.5 mm for new born babies 2 x 2.5 mm for pediatric and 2 x 4 mm for adults, so it means total 06 canula also available in a plastic tube packing.

Sturdy adjustment mount bracket.  
2 spare 2.5 V bulb.  
Pentoscope set battery handle is operating with two AA size batteries.

This Pentoscope set Packed in our new beautiful Blue Colour plastic case in same blue inner packing material.  
Spare parts also available of any item.



With 2 spare Bulbs



Mouth Mirrors Head



Otoscope Head



Sturdy Adjustment Mount Bracket



Standard illumination



Plastic ear Canulas



Two laryngeal mirrors: n. 4 and 5

Handle

Bent-arm illuminator: pin contact fitting.

Specula

**Professional Otoscope & Ophthalmology Set  
ENT Medical Diagnostic Set  
01-109**

2.5 V Standard illumination for use with 2 batteries C type.

Including:

Otoscope Head: Luir Lock fitting. Chrome plated.

Fix lens 4 x in plastic mount can be removed for introduction of instruments.

Specula sizes: 2.5, 3.5 & 4.5 mm Included,

Ophthalmoscope head: Luir Lock fitting.

Lens selection wheel permits choice of corrective lenses from 0 to - 20 and 0 to +20 diopters.

Single large circle aperture.

Bent-arm illuminator: pin contact fitting.

Nasal speculum: Chrome plated.

Two laryngeal mirrors: n. 4 and 5.

Battery handle: 28mm. For 2 cell C-batteries.

Chrome plated. For adjustment of all illumination.

Replacement bulb: 2.5 V

Diagnostic set is supplied in our new beautiful look and quality Blue Case, which protect our quality

Products from moisturizing and damaging.



With 2 spare Bulbs



Otoscope Head



Nasal speculum: Chrome plated.



Ophthalmoscope head



ENT Specula



Diagnostik-Bestestecke  
Trousses à diagnose

Juegos para diagnóstico  
Astucci diagnostici



Two laryngeal mirrors: n. 4 and 5



Bent-arm illuminator: pin contact fitting.

**Universal Otoscope & Ophthalmology Set  
ENT Medical Diagnostic Set  
01-110**

2.5 V Standard illumination for use with 2 batteries type C.

Including:

Otoscope Head: Luir Lock fitting.  
Chrome plated.

Fix lens 4 x in plastic mount can be removed for introduction of instruments.

Specula sizes: 2.5, 3.5 & 4.5 mm Included,

Ophthalmoscope head: pin contact fitting.

Lens selection wheel permits choice of corrective lenses from 0 to - 20 and 0 to +20 diopters.

Single large circle aperture.

Bent-arm illuminator: Luir Lock fitting.

Nasal speculum: Chrome plated.

Two laryngeal mirrors: n. 4 and 5.

Battery handle: 28mm. For 2 cell C-batteries.

Chrome plated. For adjustment of all illumination.

Replacement bulb: 2.5 V

The Diagnostic set is supplied in new beautiful look and quality Blue Case, which protect our quality Products from moisturizing and damaging.



Nasal speculum: Chrome plated.



ENT Specula



Otoscope Head



Ophthalmoscope Head



Handle



**Standard Otoscope + Ophthalmoscope Set**  
01-111

This Set is with 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luir / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Resuable speculas with sizes 2.5, 3.5 & 4.5 mm. A Full Functional Ophthalmoscope compact durable, light weight. Standard bulb and precision system 4 different apertures selection range of lenses -20 to +20 bright white light and brilliant image A Full Functional Ophthalmoscope Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Plastic Box Packing.



**Standard Ophthalmoscope Set**  
01-112

A Full Functional Ophthalmoscope compact durable, light weight. Standard bulb and precision system 4 different apertures selection range of lenses -20 to +20 bright white light and brilliant image Packed in Black Plastic Case.



'C' type Battery Handle



**Standard Otoscope Set**  
01-113

this otoscope is with 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luir / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Resuable speculas with sizes 2.5, 3.5 & 4.5 mm. Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Pouch Packing



Ophthalmoscope head

Diagnostik-Bestestecke  
Trousse à diagnose

Juegos para diagnóstico  
Astucci diagnostici



**Standard Otoscope Set  
01-114**

this otoscope is with 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luir / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Reusable speculas with sizes 2.5, 3.5 & 4.5 mm. Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Pouch Packing



**Otoscope + Ophthalmoscope Set  
01-115**

This otoscope is with 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luir / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Reusable speculas with sizes 2.5, 3.5 & 4.5 mm. A Full Functional Ophthalmoscope compact durable, light weight. Standard bulb and precision system 4 different apertures selection range of lenses -20 to +20 bright white light and brilliant image A Full Functional Ophthalmoscope Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Pouch Packing.




**Standard Ophthalmoscope  
01-116**

A Full Functional Ophthalmoscope. compact durable, Standard bulb and precision system 4 different apertures selection range of lenses -20 to +20 bright white light and brilliant image Handle made of Stainless Steel Material.




**Mini Ophthalmoscope  
01-117**

A Full Functional Ophthalmoscope. compact durable, light weight. Standard bulb and precision system 4 different apertures selection range of lenses -20 to +20 bright white light and brilliant image Made of plastic Material.




**Mini Ophthalmoscope  
01-118**

A Full Functional Ophthalmoscope. compact durable, light weight. Standard bulb and precision system 4 different apertures selection range of lenses -20 to +20 bright white light and brilliant image made of plastic Material.



- 01-119**  
Dissecting Set of 6 pcs.
- Dropper Oil
  - Probe Plastic Handle
  - Dressing Forceps
  - Dressing Scissors S/s
  - Hook
  - Scalpel Handle 4



- 01-120**  
Dissecting Set of 12 pcs.
- Iris Scissors
  - Scalpel Knives 4
  - Probe Single Ended
  - Scalpel Handle 3
  - Dressing Scissors S/s
  - Dressing Scissors S/b
  - Dressing Forceps
  - College Tweezer
  - Probe Straight Plastic Handle
  - Probe Angular Plastic Handle
  - Dropper Owl
  - Hook



- 01-121**  
**New Instructor's Dissecting set**  
Contain quality stainless steel Instrument Enclosed in an attractive zippered case Scissor surgical 4 1/2" scissiors fine point 4 1/2" Forceps fine point Str . Forcep fine point Cvd . Mosquito forceps Cvd 5" box lock . One piece scalpel 1 1/2" blade Surgeon handle # 4 . Mall probe and seeker 2 plastic handle str needle also 2 blade # 22 for each set.



- 01-122**  
**Dissecting Set- Advanced**  
Stainless steel, consisting of 12 instruments for accurate dissection. Supplied in zipper leather case.  
Kit consists of:  
Scissors - 2 pcs. Forceps - 2 pcs. Needle - 2 pcs.  
Artery Forceps - 1 pc. Tissue Forceps - 1 pc.  
Seeker - 1 pc. Scalpel Handle w/ Blade - 1 pc.  
Bone Cutter - 1 pc. Dental Probe & Seeker - 1 pc



- 01-123**  
**Dissecting Kit**  
Kit includes:  
1 4 1/2" Dissecting forceps with curved and serrated tips  
1 Curved dissecting needle 1 Glass pipette with rubber bulb  
1 4 1/2" Dissecting forceps straight and serrated tips  
1 Dissecting scalpel, 6" long, 1 1/2" cutting edge.  
1 4 1/2" Dissecting scissors, sharp, straight points  
1 Plastic ruler 6", 1/16" div., millimeter and centimeter



- 01-124**  
**Dissecting Set - Pro**  
Stainless steel, consisting of 20 instruments for accurate dissection. Supplied in zipper leather case.  
Kit consists of:  
Scissors - 2 pcs. Forceps - 2 pcs. Needles - 4 pcs. Scalpel - 1 pc.  
B.P. Handle w/ Blade - 1 pc. Razor - 1 pc. Ruler 6" - 1 pc.  
Magnifying Glass - 1 pc. Glass Slide - 1 pc. Chain and Hook - 1 pc.  
Cover Slips - 1 pc. Blow Pipe - 1 pc. Needle Wooden Handle - 1 pc.  
Dropper Glass - 1 pc. Brush - 1 pc.



- 01-125**  
**Dissecting Set - Popular**  
Stainless steel, consisting of 14 instruments for accurate dissection. Complete in zipper leather case.  
Kit consists of:  
Scissors - 2 pcs. Forceps - 2 pcs. Needles - 4 pcs. Scalpel - 1 pc.  
B.P. Handle w/ Blade - 1 pc. Razor - 1 pc. Magnifying Glass - 1 pc.  
Chain and Hook - 1 pc. Blow Pipe - 1 pc.



- 01-126**  
**Dissecting Set - Senior**  
Stainless steel, consisting of 9 instruments for accurate dissection. Supplied in zipper leather case.  
Kit consists of:  
Scissors - 2 pcs. Forceps - 2 pcs. Needle - 2 pcs. Scalpel - 1 pc.  
Scalpel Handle w/ Blade - 1 pc. Ruler 6" - 1 pc.

Diagnostik-bestecke  
Estuches de diagnostico

Diagnostic Setts  
Combinazioni per diagnostica



**01-127**  
Dissecting Set



**01-128**  
Dissecting Set



**01-129**  
Dissecting Set



**01-130**  
Dissecting Set



**01-131**  
Dissecting Set



**01-132**  
Dissecting Set



**Jameson Caliper**  
01-133 8 cm



**Townley Caliper**  
01-134 10 cm



**Vernier Caliper**  
01-135 15 cm



**Vernier Caliper**  
01-136 15 cm



**Carroll Pocket Goniometer**  
01-137 14.5cm



**Möeltgen Flexometer**  
01-138 24cm



**Bertillon Measuring Device**  
01-139 32 cm



**Castroviejo-Epker Caliper**  
01-140 18 cm  
Measuring range 40 mm, straight



**Castroviejo-Epker Caliper**  
01-141 18 cm Angled  
Measuring range 40 mm



**Castroviejo Bone Caliper**  
01-142 18cm  
Straight



**Castroviejo Bone Caliper Curved**  
01-143 18cm



**Castroviejo-Epker Caliper**  
01-144 9 cm Str.  
Measuring range, 40 mm  
01-145 9 cm Cvd.  
Measuring range, 40 mm



**Castroviejo Eye View Caliper**  
01-146 8 cm, Str.  
Measuring range, 20 mm  
01-147 8 cm, Cvd.  
Measuring range, 20 mm



**Castroviejo Eye View Caliper**  
01-148 8 cm  
Measuring range 20 mm



**Tonometer**  
01-149 Set

